

## PRODUCT DATASHEET

Model: **PDRNGDEVO43**

Description: **Fuji D-EVO 43x43 Mobile DR Panel Holder**

Revision: May 2017

### PRODUCT ATTRIBUTES:

- Holds DR Panel during exams.
- Provides handle and lock mechanism to protect panel during handling.
- Fits Fuji D-EVO Advanced 43x43 DR Panel.

**Reina**Imaging



### MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	HDP	98
Overlay	Polycarbonate Film	2

### DIMENSIONS AND PERFORMANCE

Product Dimension: 20.86 x 18.81x 1.096 in. (52.98 x 47.77 x 2.78 cm)

Weight: 1.95 lbs. (0.8845 kg)

Weight Tolerance: 400 lbs. (181.4 kg)

*Not designed for weight bearing exams. Not to be used without an image receptor.*

### FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

### STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

### REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

6107 Lou Ave.

Crystal Lake, IL 60014 USA

### AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

### DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | [reinaimaging.com](http://reinaimaging.com)