

Product Datasheet

Model: PDRNG3543EZ
Description: MOBILE DR PANEL HOLDER FOR THALES 3543EZ
Revision: SEPTEMBER 2013

Product Attributes:

- Holds DR Panel during exams.
- Protects DR Panel during exams and storage.
- Provides additional protection for DR Panel.

Dimensions

Product Dimension: 18.9 x 18.8 x 1in. (48.1 x 47.8 x 2.54cm)
Weight: 1.7 lb. (0.77kg)

Material Composition

Part	Material	Approximate %
Shell	HDP	98
Overlay	Lexan Film	2

Weight Tolerance

Note: This was not designed for weight-bearing studies and should never be used without the DR Panel installed.

Weight Tolerance: 400lb. (181.4kg)

Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.
UL Rating: UL94V-2.

Biocompatibility

This model is a biocompatible device per ISO 10993.

Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

Registrations and Compliance

All models listed above are made in the USA and registered to:

Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA

Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1