

PRODUCT DATASHEET

Model: PDRNG810C

Description: CXDI-810C MOBILE DR PANEL HOLDER

Revision: March 2021



Product Attributes

- Encases grid to protect against damage
- Holds DR Panel during exams
- Provides additional protection for DR Panel

Material Composition

Material	%
Shell, HDP	98%
Overlay, Polycarbonate Film	2%

DIMENSIONS

Product Dimension: 15x16x1in (38.1x40.64x2.54cm)

Weight: 1.55lb. (.67kg)

Weight Tolerance: 400 lbs. (181.4 kg)

Not designed for weight bearing exams. Not to be used without an image receptor.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA